

Summary Report – Run 112 p63/AMACR/HMWCK

Assessors: E Belanger, J Garratt, B Gilks, J Won

Overview

Expression of AMACR combined with loss of expression of basal cell markers, HMWCK and p63, enhance diagnostic accuracy of prostate cancer. The survey consisted of eight (8) 2mm tissue cores: five (5) prostate adenocarcinomas and three (3) benign prostate cores. Overall assessment of immunohistochemical staining in all eight cores combined was used to evaluate staining of each marker. The adenocarcinoma in Core 7 was cut through for some participants. It was also noted that Cores 2 and 7 may not have been as well-fixed as other cores in the tissue microarray, leading to generally weaker staining for p63 in particular. The expected staining criteria are described below:

p63

- Nuclear staining in benign basal cell nuclei of moderate intensity
- No nuclear staining in carcinoma
- Sufficient staining intensity
- Minimal to no non-specific staining

AMACR (α -methylacyl-coenzyme A racemase)

- Granular cytoplasmic staining of tumour with luminal accentuation in glands of adenocarcinoma that lack p63
- At most, only weak focal staining in benign
- Sufficient staining intensity

Prostate basal epithelial keratin

- Nuclear staining in benign basal cell nuclei of moderate intensity
- No nuclear staining in carcinoma
- Sufficient staining intensity
- Minimal to no non-specific staining

Results

The technical quality of staining for participants was excellent. Participant-specific feedback is below:


Lab ID	p63		AMACR		HMWCK	
	IHC Status*	Comment	IHC Status*	Comment	IHC Status*	Comment
101	Optimal		Optimal		Optimal	
102	Optimal		Optimal		Optimal	
103	Optimal		Optimal		Optimal	
111	Optimal	Slightly weak	Optimal		Optimal	
112	Optimal		Optimal	Slightly more staining of benign prostate than other labs	Optimal	
113	Optimal		Optimal	Stromal background	Optimal	
120	Optimal		Optimal		Optimal	
132	Optimal		Optimal		Optimal	
147	Optimal		Optimal	Stromal background	Optimal	
149	Optimal		Optimal		Optimal	
151	Adequate	Cores 2 and 7 poor staining due to weak staining	Optimal		Optimal	
183	Optimal		Optimal		Optimal	
198	Optimal		Optimal	Slight stromal background	Optimal	
207	Optimal	Slightly weak	Optimal		Optimal	
230	Optimal		Optimal	Nice staining	Optimal	
231	Optimal		Optimal		Optimal	
234	Optimal		Optimal	Slightly weak	Optimal	

*based on CPQA assessor consensus

Garrattogram after CPQA assessment (provided as an overall assessment of all tissue microarray cores):

Lab ID	102	103	111	112	113	120	132	147	149	151	183	198	207	230	231	234
p63	○	○	○	○	○	○	○	○	○	A	○	○	○	○	○	○
AMACR	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
HMW CK	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○

Lab 102 (Optimal)



Lab 151 (Adequate; weak staining)




Figure 1. Representative images of the qualitative variability of p63 staining.

Lab 230
(Optimal)

Lab 112
(Optimal)

Lab 113
(Optimal; stromal staining)

Lab 234
(Optimal; slightly weak)




Figure 2. Representative images of the qualitative variability of AMACR staining in benign prostate and prostate adenocarcinoma.

Supplementary Tables 1-3 summarize the reported staining protocols for p63, AMACR and prostate basal epithelial keratins, respectively, which can be referred to during validation or optimization of a staining protocol. Supplementary Table 4 provides the definitions of IHC Status and recommended participant action. Your regular participation in CPQA is greatly appreciated and we look forward to continuing to work with you and the Canadian Association of Pathologists – Association Canadienne des Pathologistes.

This report has been updated with scanned images that were acquired using a NanoZoomer SQ that has been graciously loaned to the CPQA-AQCP by Quorum Technologies and Hamamatsu.

Table S1. Reported p63 staining protocols.

Lab ID	Ag Retrieval Method	Time for Ag Retrieval (min)	Ab Clone	Ab Dilution	Ab Supplier/ Vendor	Ab Lot No.	Time for Ab Incubation (min)	Detection System	Amplification (Y/N)	Enhancement (Y/N)	Chromogen
102	DAKO PT - HIGH PH	20"	4A4	1:50	BIOCARE	52219	30" RT	DAKO ENVISION FLEX	NO	YES CUSO4	DAB+
103	CC1	36	4A4	Pre	Ventana	F22650	32	Ultra view	Y	Y copper	DAB
112	Bond ER2 pH 9.0	30 minutes	4A4	RTU	BioCare Medical	80819	60 minutes	ChromoPlex 1 Dual Detection	no	no	DAB
113	High pH	30min	DAK-p63	1/100	Dako	2007039	27.5min	Dako Envision Flex	10min Mouser Linker	N	DAB
120	HIER Waterbath	20	DAK-p63	RTU	Dako	20066310	20	Envision Flex+	N	N	DAB
132	High pH	20	DAK-P63	RTU	Dako	20066311	20	Envision flex	N	N	DAB
147	ER PH 9	20	13H4	50	ZETA CORPORAT ION	1708B/1902 C-1	15	LEICA REFINE KIT	N	N	DAB
149	high pH OMNIS	20 min at 97 C	DAK-p63	RTU	Dako Agilent	20066314	20	EnVision Flex OMNIS	No	No	DAB
151	BUFFER 9.0	20 MIN	4A4	1:200	ZETA	181221	15 MIN	BOND REFINE	N	N	DAB
183	ER2	30	CK5/14+ P63+P50 4S	RTU	BIOCARE	21519	15	CHROMOPLEX	N	N	DAB
198	HIER high pH	30 min @ 97deg	DAK-p63	RTU	Agilent	20073135	15 min	Envision Flex	N	N	DAB
207	on line	64	4A4	prediluted	Ventana	F22650	32	ULTRAVIEW	N	Y	DAB
230	HIER	40	A4A	predilute	Roche	F18327	32	Optiview	N	N	DAB
231	HIER	56 mins	4A4	RTU	ROCHE/VE NTANA	F00224	24 mins	OPTIVIEW (VENTANA)	N	Y	DAB
234	Omnis on board	30	DAK-p63	RTU	Agilent	20073135	10	Envision	Rabbitinker	n	DAB

Table S2. Reported AMACR staining protocols.

Lab ID	Ag Retrieval Method	Time for Ag Retrieval (min)	Ab Clone	Ab Dilution	Ab Supplier/ Vendor	Ab Lot No.	Time for Ab Incubation (min)	Detection System	Amplification (Y/N)	Enhancement (Y/N)	Chromogen
102	DAKO PT - HIGH PH	20"	13H2	1:80	DAKO	10148602	30" RT	DAKO ENVISION FLEX	NO	YES CUSO4	DAB+
103	CC1	36	13H4	1/100	Dako	10143476	24	Ultra View	No	Y copper	Dab
112	Bond ER 2 pH 9.0	30 minutes	N/A (polyclonal)	RTU	BioCare Medical	80819	60 minutes	ChromoPlex 1 Dual Detection	no	no	Fast Red
113	High pH	30min	13H4	Pre-dilute	Dako	11102752	20min	Dako Envision Flex	10min Rabbit Linker	N	DAB
120	HIER Waterbath	20	13H4	RTU	Dako	10147735	20	Envision Flex+	N	N	DAB
132	High pH	20	13H4	1:50	Dako	10136945	30	Envision flex	N	N	DAB
147	ER PH 9	20	13H4	50	ZETA CORPORAT ION	1708B/1902 C-1	15	LEICA REFINE KIT	N	N	DAB
149	high pH OMNIS	20 min at 97 C	13H4	RTU	Dako Agilent	10146449	15	EnVision Flex OMNIS	Yes	No	DAB
151	BUFFER 9.0	20min	13H4 +4A4	1:150	ZETA	1703A/1072 A-1	15 MIN	BOND REFINE	N	N	DAB
183	ER2	30	CK5/14+ P63+P50 4S	RTU	BIOCARE	21519	15	CHROMOPLEX	N	N	REFINE RED
198	HIER High pH	30 min 97deg	13H4	RTU	Agilent	11102752	20	Envision Flex	Y	N	Magenta Envision
207	ON LINE	64	P504S	1/200	BIOCARE	72318	44	ULTRAVIEW	N	Y	DAB
230	HIER	40	P504S +4A4	predilute	Biocare	VP201G	32	Optiview	N	N	DAB
231	HIER	52 mins	polyclonal	1/100	intermedico (biocare)	70319	32 mins	ULTRAVIEW (VENTANA)	N	Y	DAB
234	Omnis on board	30	13H4	5000	Agilent	10136945	15	Envision	Rabbitinker	n	DAB

Table S3. Reported prostate basal epithelial keratins staining protocols.

Lab ID	Ag Retrieval Method	Time for Ag Retrieval (min)	Ab Clone	Ab Dilution	Ab Supplier/ Vendor	Ab Lot No.	Time for Ab Incubation (min)	Detection System	Amplification (Y/N)	Enhancement (Y/N)	Chromogen
102	DAKO PT - HIGH PH	20"	D5/16 B4	Predilute	DAKO	20067284	30" RT	DAKO ENVISION FLEX+	YES	YES CUSO4	DAB+
103	CC1	36	34BE12	Pre	Ventana	F19791	16	Ultra view	N	Y copper	DAB
112	Bond ER2 pH 9.0	30 minutes	XM26 & LLOO2	RTU	BioCare Medical	80819	60 minutes	ChromoPlex 1 Dual Detection	no	no	DAB
113	High pH	30min	XM26	1/200	Leica	6065470	27.5min	Dako Envision Flex	10min Mouse Linker	N	DAB
120	HIER Waterbath	20	34BE12	RTU	Dako	10146592	20	Envision Flex+	N	N	DAB
132	High pH	20	34BE12	RTU	Dako	10149455	20	Envision flex	N	N	DAB
147	ER PH 9	20	34BE12	100	DAKO	10147751	15	LEICA REFINE KIT	N	N	DAB
149	high pH OMNIS	20 min at 97 C	34BE12	RTU	Dako Agilent	10138647	20	EnVision Flex OMNIS	Yes	No	DAB
151	BUFFER 9.0	20 MIN	34BE12	1:200	DAKO	10147751	20 MIN	BOND REFINE	N	N	DAB
183	ER2	30	CK5/14+ P63+P50 4S	RTU	BIOCARE	21519	15	CHROMOPLEX	N	N	DAB
198	HIER high pH	30 min @97 deg	34BE12	RTU	Agilent	11086028	25 min	Envision Flex	Y	N	DAB
207	ON LINE	64	34BE12	PREDILUT ED	VENTANA	F13474	16	ULTRAVIEW	N	Y	DAB
230	HIER	40	34BETA E12	predilute	Roche	F04650	32	Optiview	N	N	DAB
231	HIER	36 mins	34BE12	RTU	ROCHE/VENTANA	F19791	32 mins	ULTRAVIEW (VENTANA)	N	Y	DAB
234	Omnis on board	30	34BE12	RTU	Agilent	10142752	25	Envision	Mouselinker	n	DAB

Table S4. Proficiency Testing Definitions of IHC Status.

IHC Status	Definition	Proficiency Testing Performance
Optimal	All expected targets are identified appropriately and demonstrate the expected staining intensity. Absence of non-specific staining (no background staining).	PASS
Adequate	All targets are identified, but intensity of staining is weaker than optimal or there is false-positive staining which does not interfere with interpretation.	PASS
Sub-optimal	None or only some targets are identified OR all targets are identified, but false-positive staining may interfere with interpretation.	PASS, Conditionally¹
Failed	The staining was considered to be of such poor quality that accurate readout of the test is unlikely or impossible.	FAIL²
Unsatisfactory	Technical issue (e.g. unsuitable antibody selection, etc.)	N/A

¹ - A one-time sub-optimal performance qualifies for a "Pass" result. Two successive "sub-optimal" results will be designated as a "Fail".

^{1,2} - Please contact the CPQA for assistance and, if necessary, inform your regional regulatory body as per the terms of your laboratory's accreditation provider.